

Product Data Sheet

Cat# SP-101110-5

Description: [Ala5, β-Ala8]-Neurokinin A (4-10)[Ala5, β-Ala8]-Neurokinin A (4-10)

Sequence: Asp-Ala-Phe-Val-β-Ala-Leu-Met-NH₂; MW: 765

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-101110-5 101110A