

Product Data Sheet

Cat# **SP-101113-1**

Description: [D-Arg25]-Neuropeptide Y, human, rat; [D-Arg25]-NPY, human, rat; [D-Arg25]-Neuropeptide Y, human, rat; [D-Arg25]-NPY, human, rat

Sequence: Tyr-Pro-Ser-Lys-Pro-Asp-Asn-Pro-Gly-Glu-Asp-Ala-Pro-Ala-Glu-Asp-Met-Ala-Arg-Tyr-Tyr-Ser-Ala-Leu-D-Arg-His-Tyr-Ile-Asn-Leu-Ile-Thr-Arg-Gln-Arg-Tyr-NH₂; MW: 4271.80

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-101113-1 101110A