

---

<b>Cat#</b>	<b>SP-101264-1</b>
<b>Description:</b>	[Des-Gln16]-PACAP (6-27), amide, human, ovine, rat[Des-Gln16]-PACAP (6-27), amide, human, ovine, rat
<b>Sequence:</b>	Phe-Thr-Asp-Ser-Tyr-Ser-Arg-Tyr-Arg-Lys-Met-Ala-Val-Lys-Lys-Tyr-Leu-Ala-Ala-Val-Leu-NH <sub>2</sub> ; MW: 2510
<b>Size:</b>	1 mg
<b>Purity:</b>	>95%
<b>Store:</b>	Desiccated at -20oC.

---

*All peptides are for in vitro research use only.*

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.  
[http://www.4adi.com/commerce/info/showpage.jsp?page\\_id=1088&category\\_id=2427](http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427)

**SP-101264-1**                      101110A