

Product Data Sheet

Cat# SP-101512-1

Description: [Ser25] - PKC (19 - 31), biotinylated[Ser25] - PKC (19 - 31), biotinylated

Sequence: Lys(Biotin)-Arg-Phe-Ala-Arg-Lys-Gly-Ser-Leu-Arg-Gln-Lys-Asn-Val; MW 1914.3

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-101512-1 101110A