

Product Data Sheet

Cat# **SP-101680-1**

Description: Ac-Angiotensinogen (1-14), porcineAc-Angiotensinogen (1-14), porcine

Sequence: Ac-Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu-Leu-Val-Tyr-Ser (MW: 1801.1)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-101680-1 101110A