

Product Data Sheet

Cat# SP-101851-5

Description: [Tyr1] Adipokinetic Hormone, locust[Tyr1] Adipokinetic Hormone, locust

Sequence: Tyr-Leu-Asn-Phe-Thr-Pro-Asn-Trp-Gly-Thr-NH₂; MW 1211.5

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-101851-5 101110A