

Product Data Sheet

Cat# **SP-102106-5**

Description: [Ala9] Autocamtide 2; Autocamtide-2-Related Inhibitory Peptide[Ala9] Autocamtide 2; Autocamtide-2-Related Inhibitory Peptide

Sequence: Lys-Lys-Ala-Leu-Arg-Arg-Gln-Glu-Ala-Val-Asp-Ala-Leu; MW: 1497.7

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-102106-5 101110A