

0

Cat# SP-51029-1

Description: MAGE-3 Antigen (271-279) peptide

Sequence: H-Phe-Leu-Trp-Gly-Pro-Arg-Ala-Leu-Val-OH; MW 158.3

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.



SP-51029-1 0

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com