

Product Data Sheet

Cat# SP-53014-5

Description: Human CMV pp65 (495-503) (AA: Asn-Leu-Val-Pro-Met-Val-Ala-Thr-Val) (MW: 943.18)

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-53014-5 101115A