

Product Data Sheet

Cat# **SP-53024-5**

Description: Melanocyte Protein PMEL 17 (185-193) (human, bovine, mouse) (AA: Ile-Thr-Asp-Gln-Val-Pro-Phe-Ser-Val) (MW: 1005.14)

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-53024-5 101115A