
Product Data Sheet

Cat# **SP-54677-1**

Description: Histatin 3 (H3) (AA: Asp-Ser-His-Ala-Lys-Arg-His-His-Gly-Tyr-Lys-Arg-Lys-Phe-His-Glu-Lys-His-His-Ser-His-Arg-Gly-Tyr-Arg-Ser-Asn-Tyr-Leu-Tyr-Asp-Asn) (MW: 4062.44)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

SP-54677-1 120703A