

Product Data Sheet

Cat# **SP-54841-1**

Description: LEAP-2 (38 - 77) (Human) (AA: Met-Thr-Pro-Phe-Trp-Arg-Gly-Val-Ser-Leu-Arg-Pro-Ile-Gly-Ala-Ser-Cys-Arg-Asp-Asp-Ser-Glu-Cys-Ile-Thr-Arg-Leu-Cys-Arg-Lys-Arg-Arg-Cys-Ser-Leu-Ser-Val-Ala-Gln-Glu) (MW: 4585.36)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-54841-1 101115A