

Cat# SP-55196-1 Plasminogen Activator Inhibitor 1 (PAI-1) hexapeptide

Size: 0.5 mg

Names	Plasminogen Activator Inhibitor 1 / PAI-1, hexapeptide (human)
Sequence	Ala - Arg - Met - Ala - Pro - Glu
MW	673.33
Purity	>95%
Solubility	Soluble in water
Notes	For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. Do not attempt to weigh out smaller portions

References Blair-West JR (1971) J Clin Endocrinol Metab. 32:575-8

SP-55196-1

rev 90921A