

Product Data Sheet

Cat# **SP-55278-05**

Description: Atrial Natriuretic Peptide (1-28), rat (AA: Ser-Leu-Arg-Arg-Ser-Ser-Cys-Phe-Gly-Gly-Arg-Ile-Asp-Arg-Ile-Gly-Ala-Gln-Ser-Gly-Leu-Gly-Cys-Asn-Ser-Phe-Arg-Tyr
-Tyr
□Disulfide bridge:Cys7-Cys23□(MW: 3062.47)

Size: 0.5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-55278-05 101115A