

Product Data Sheet

---

**Cat#**            **SP-86706-1**

**Description:** Histone H3 (116–136), C116–136 (AA: Lys-Arg-Val-Thr-Ile-Met-Pro-Lys-Asp-Ile-Gln-Leu-Ala-Arg-Arg-Ile-Arg-Gly-Glu-Arg-Ala) (MW: 2508.06)

**Size:**            1 mg

**Purity:**         >95%

**Store:**           Desiccated at -20oC.

---

SP-86706-1            120703A