

Cat# **SP-89550-1**

Description: [Tyr⁰]-C-Type Natriuretic Peptide (32-53) (human, porcine, rat) [Tyr-Gly-Leu-Ser-Lys-Gly-Cys-Phe-Gly-Leu-Lys-Leu-Asp-Arg-Ile-Gly-Ser-Met-Ser-Gly-Leu-Gly-Cys (Disulfide bridge: Cys7-Cys23); MW 2360.82]

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

SP-89550-1 120703A