

Product Data Sheet

Cat # SPL17-R

Recombinant (E. coli) Syphilis antigen (17 kDa), purified

Size: 100 ug

Syphilis is a sexually transmitted disease caused by the spirochetal bacterium *Treponema pallidum* subspecies *pallidum*. The route of transmission of syphilis is almost always through sexual contact, although there are examples of congenital syphilis via transmission from mother to child in utero or at birth. Syphilis is a disease caused by Spirochete bacterium called *Treponema pallidum* (TP). If untreated, the organisms move throughout the body and can cause damage to many organs, making syphilis a life-threatening disease if not treated early enough. People who have been infected with Syphilis experience different symptoms during the 3 stages of the disease. Early, which is defined by the presence of the chancre at the site of inoculation. Syphilis may be further divided into primary, secondary, and early latent syphilis; late syphilis includes late latent and the various forms of tertiary Syphilis. The serological response to syphilis involves production of antibodies to a wide range of antigens, including non-specific antibodies and specific anti-TP antibodies. The first detectable response to infection is the production of specific antitreponemal IgM, which can be detected within 4 to 7 days after the chancre appears and until the end of the second week of infection; antitreponemal IgG appears at about four weeks later. By the time that symptoms develop, most patients have detectable IgG and IgM.

important antigens and potential virulence factors have been identified in *T. pallidum*. One such antigen, with a molecular mass of 47-kilodaltons (kDa), has been shown to be surface-associated, abundant, and highly immunogenic in both human and experimental syphilis. Monoclonal antibodies to the 47-kDa antigen possess complement-dependent treponemicidal activity in both the *T. pallidum* immobilization test and the in vitro-in vivo neutralization test.

**Formulation:**

Syphilis 17 kDa protein is produced in *E. coli* and purified (>95%). The protein is supplied in 50mM Tris, (pH7.2) with 0.05% NaN<sub>3</sub> in liquid or lyophilized in the same buffer (see lot specific cocn on the vial). It is recommended to reconstitute the lyophilized protein in sterile water at not less than 100µg/ml, which can then be further diluted to other aqueous solutions. Lyophilized protein although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution, protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

**References:** Baker-Zander SA (1985) *J. Inf. Immun. Dis.* 151, 264-272; Hanff PA (1983) *Inf. Immun.* 40, 825-828; . 129, 1287-1291; Neal R (1988) *Inf. Immun.* 56, 71-78

Related items

Catalog#	ProdDescription	Quantity/SIZE
4910	Anti-Tp ( <i>Treponema pallidum</i> ) Syphilis ELISA kit, Semi-Quantitative	1 kit
4915-R	Anti-Tp ( <i>Treponema pallidum</i> ) Syphilis rapid tests strips	100 pk
SPL15-R	Recombinant ( <i>E. coli</i> ) Syphilis antigen (15kDa), purified	100 ug
SPL17-R	Recombinant ( <i>E. coli</i> ) Syphilis (17kDa), purified	100 ug
SPL15-R	Recombinant ( <i>E. coli</i> ) Syphilis (47kDa), purified	100 ug
SPL17-R		100809A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)