

ELISA Reagent Set for Melon necrotic spot virus (MNSV)



Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The MNSV ELISA is an effective tool for the detection of *Melon necrotic spot virus* in cucurbit plant tissues.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)
- [Tech Note](#)

Pack Sizes

Reagent Set for Melon necrotic spot virus (MNSV), 96 testwells SRA 12402/0096

Reagent Set for Melon necrotic spot virus (MNSV), 500 testwells SRA 12402/0500

Reagent Set for Melon necrotic spot virus (MNSV), 1000 testwells SRA 12402/1000

Reagent Set for Melon necrotic spot virus (MNSV), 5000 testwells SRA 12402/5000