

ELISA Reagent Set for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph)

This ELISA test is a qualitative serological assay for the detection of *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), the causal agent of common blight of bean in bacterial cultures, bean leaves, and bean seeds. Xaph is a member of the *Xanthomonas* genus known for its many species and pathovars infecting over 400 plant species.

Xaph was formerly known as *Xanthomonas campestris* pv. *phaseoli*

Included:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)

Option

Reagent Set for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 96 testwells

SRA 36101/0096

Reagent Set for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 500 testwells

SRA 36101/0500

Reagent Set for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 1000 testwells

SRA 36101/1000

Reagent Set for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 5000 testwells

SRA 36101/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com