

ELISA Reagent Set for *Plantago asiatica* mosaic virus (PIAMV)



This ELISA test is a qualitative serological assay for the detection of *Plantago asiatica* mosaic virus (PIAMV) in ornamental bulb crop leaves. PIAMV is a member of the Potexvirus genus known for their non-enveloped, flexuous, filamentous-shaped virus particles.

If you would prefer to add a field-deployable assay to your testing toolkit, please check out Agdia's [PIAMV ImmunoStrips®](#)

<input type="checkbox"/>	<p><i>included:</i></p> <ul style="list-style-type: none"> • Capture antibody • Alkaline phosphatase enzyme conjugate • 96-well microtiter plates, solid* (*strip in single plate set) • User Guide <p>Option</p>
<input type="checkbox"/>	<p>Reagent Set for <i>Plantago asiaticus</i> mosaic virus (PIAMV), 96 testwells SRA 91500/0096</p>
<input type="checkbox"/>	<p>Reagent Set for <i>Plantago asiaticus</i> mosaic virus (PIAMV), 500 testwells SRA 91500/0500</p>
<input type="checkbox"/>	<p>Reagent Set for <i>Plantago asiaticus</i> mosaic virus (PIAMV), 1000 testwells SRA 91500/1000</p>
<input type="checkbox"/>	<p>Reagent Set for <i>Plantago asiaticus</i> mosaic virus (PIAMV), 5000 testwells SRA 91500/5000</p>

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com