

ImmunoStrip® for mBt-Cry3A

Part Number: STX 06700/0050

Test Label: NA

Test Format: ImmunoStrip®

Capture Reagent: Monoclonal

Detection Reagent: Monoclonal

Intended Use:

Especially designed for plant breeders, the Agdia mBt-Cry3A ImmunoStrip® effectively detects the mBt-Cry3A protein from Agrisure® RW (MIR6042) hybrid or inbred seed corn and leaves.

Once sample is prepared, test completion times are usually between 3 to 5 minutes.

The mBt-Cry3A ImmunoStrip® is applicable for use in the lab or the field, but special emphasis was placed on making a test useful for plant breeders under field conditions. All Agdia ImmunoStrips® are rugged, reliable and ready for field use.

Includes:

- ImmunoStrips®
- SEB4 sample buffer not included, see related products
- [User Guide](#)

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com