

## ImmunoStrip® for Potato virus Y (PVY)



**Test Label:** NA  
**Test Format:** ImmunoStrip®  
**Capture Reagent:** Polyclonal  
**Detection Reagent:** Polyclonal

### Intended Use:

The PVY ImmunoStrip® detects of all strains of PVY in foliage and sprouts of potato.

The test is highly sensitive (comparable to ELISA methods) and can be completed in as little as a few minutes with no equipment or expertise.

### Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB4 buffer
- [User Guide](#)

### ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB4 extraction buffer purchased separately

### Option

ImmunoStrip® for Potato virus Y (PVY), 5 strips and buffer filled bags  
ISK 20001/0005

ImmunoStrip® for Potato virus Y (PVY), 25 strips and buffer filled bags  
ISK 20001/0025

ImmunoStrips® only for Potato virus Y (PVY) \*\*\*BUFFER FILLED BAGS  
NOT INCLUDED\*\*\*  
STX 20001/0025

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)