

## ImmunoStrip® for Pepper mild mottle virus (PMMoV)



The *Pepper mild mottle virus* (PMMoV) ImmunoStrip® is used to detect the presence of PMMoV in pepper leaves and seeds. This ImmunoStrips® has been shown to reliably detect 1 PMMoV-infected seed in a subsample of 249 PMMoV-free pepper seeds. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®

### Option

ImmunoStrip® for Pepper mild mottle virus (PMMoV), 5 strips and buffer filled bags  
ISK 37502/0005

ImmunoStrip® for Pepper mild mottle virus (PMMoV), 25 strips and buffer filled bags  
ISK 37502/0025

ImmunoStrip® for Pepper mild mottle virus (PMMoV), 25 strips ONLY \*\*\* BUFFER FILLED BAGS NOT INCLUDED \*\*\* STX 37502/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)