

## ImmunoStrip® for Potato virus S (PVS)



**Test Label:** NA  
**Test Format:** ImmunoStrip  
**Capture Reagent:** Polyclonal  
**Detection Reagent:** Polyclonal

### Intended Use:

The PVS ImmunoStrip detects all strains of PVS in potato leaf tissue.

The test is highly sensitive (comparable to ELISA methods) and can be completed in as little as a few minutes with no equipment or expertise making them perfect for use in the field or greenhouse.

### Kit (ISK) Includes:

- ImmunoStrips
- Sample bags containing SEB4 buffer
- [User Guide](#)

### ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB4 extraction buffer purchased separately

Option
ImmunoStrip® for Potato virus S (PVS), 5 strips and buffer filled bags ISK 40000/0005
ImmunoStrip® for Potato virus S (PVS), 25 strips and buffer filled bags ISK 40000/0025
ImmunoStrips® only for Potato virus S (PVS) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 40000/0025

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)