

ImmunoStrip® for Tobacco mosaic virus (TMV)

The *Tobacco mosaic virus* (TMV) ImmunoStrip® is used to detect the presence of TMV in many different crops including vegetables, fruits, and ornamentals. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

The TMV ImmunoStrip® cross-reacts with many other Tobamoviruses. For more technical information, [click here](#).

Tobacco mosaic virus is an extremely stable virus and can easily infect adjacent plants and contaminate surfaces such as benches, pots, tools, etc.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for Tobacco mosaic virus (TMV), 5 strips & buffer bags ISK 57400/0005
<input type="checkbox"/>	ImmunoStrip® for Tobacco mosaic virus (TMV), 25 strips & buffer bags ISK 57400/0025
<input type="checkbox"/>	ImmunoStrips® only for Tobacco mosaic virus (TMV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 57400/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com