

ImmunoStrip® for Citrus tristeza virus (CTV)



The *Citrus tristeza virus* (CTV) ImmunoStrip® is used to detect the presence of CTV in citrus leaves, young shoots, and petioles. The CTV ImmunoStrips® are **not** recommended for woody tissue. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab. This test has been proven to detect over 87 CTV isolates from all over the world, however, it cannot differentiate mild and severe strains.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

Option

ImmunoStrip® for Citrus tristeza virus (CTV), 5 strips & buffer filled bags
ISK 78900/0005

ImmunoStrip® for Citrus tristeza virus (CTV), 25 strips & buffer filled bags
ISK 78900/0025

ImmunoStrip® only for Citrus tristeza virus (CTV), ***BUFFER FILLED BAGS NOT INCLUDED***
STX 78900/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com