

ImmunoStrip® for Calibrachoa mottle virus (CbMV)



Test Label: N/A
Test Format: ImmunoStrip®
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

The CbMV ImmunoStrip® is a rapid detection tool for the presence of *Calibrachoa mottle virus* in plant foliage.

ImmunoStrip® tests require no special equipment or tools to run, making them perfect for use anywhere. Results are obtained in as little as a few minutes.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option

ImmunoStrip® for Calibrachoa mottle virus (CbMV), 5 strips and buffer filled bags
ISK 83000/0005

ImmunoStrip® for Calibrachoa mottle virus (CbMV), 25 strips and buffer filled bags
ISK 83000/0025

ImmunoStrips® only for Calibrachoa mottle virus (CbMV) ***BUFFER FILLED BAGS NOT INCLUDED***
STX 83000/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com