

Agdia ImmunoStrip® for Alfalfa mosaic virus (AMV)



The *Alfalfa mosaic virus* (AMV) ImmunoStrip® is used to detect the presence of AMV in ornamental and vegetable crops. The AMV ImmunoStrip® has been designed to detect all strains and isolates of AMV and has experimentally detected both Subgroup I and II. AMV is a member of the Alfamovirus genus known for their bacilliform, with a constant diameter of 18 nm, viral particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for AMV, 5 strips and buffer filled bags ISK 87601/0005
<input type="checkbox"/>	ImmunoStrip® for AMV, 25 strips and buffer filled bags ISK 87601/0025
<input type="checkbox"/>	ImmunoStrip® for AMV ***BUFFER FILLED BAGS NOT INCLUDED*** STX 87601/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com