

ELISA Reagent Set for Hibiscus chlorotic ringspot virus (HCRSV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The *Hibiscus chlorotic ringspot virus* ELISA is used to detect the presence of HCRSV in Hibiscus.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)

Reagent Set for HCRSV, 96 testwells, SRA 12300/0096

Reagent Set for HCRSV, 500 testwells, SRA 12300/0500

Reagent Set for HCRSV, 1000 testwells, SRA 12300/1000

Reagent Set for HCRSV, 5000 testwells, SRA 12300/5000



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com