

## ELISA Reagent Set for High plains virus (HPV)

**Test Label:** Alkaline phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

This test is appropriate for testing leaves and seed of corn and leaves of wheat and barley for *High plains virus* (HPV).

This test shows no cross-reactivity with *Wheat streak mosaic virus* which is often found as a co-infection with HPV. This test does not react with the newly emerging HPV variants.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for High plains virus (HPV), 96 testwells,SRA 17200/0096

Reagent Set for High plains virus (HPV), 500 testwells,SRA 17200/0500

Reagent Set for High plains virus (HPV), 1000 testwells,SRA 17200/1000

Reagent Set for High plains virus (HPV), 5000 testwells,SRA 17200/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)