

ELISA Reagent Set for Potato virus Y (PVY)

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This test is used for detection of all strains of *Potato virus Y* (PVY) in foliage and sprouts of potato. It can also detect the virus in several other solanaceous crops.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Potato virus Y (PVY), 96 testwells, SRA 20001/0096

Reagent Set for Potato virus Y (PVY), 500 testwells, SRA 20001/0500

Reagent Set for Potato virus Y (PVY), 1000 testwells, SRA 20001/1000

Reagent Set for Potato virus Y (PVY), 5000 testwells, SRA 20001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com