

ELISA Reagent Set for Barley yellow dwarf virus-pav (BYDV-PAV)

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal with Anti-animal conjugate

Intended Use:

This ELISA is a qualitative serological assay intended for detection of *Barley yellow dwarf virus-PAV* (BYDV-PAV) in plant foliage. The test is best suited for testing leaf and root tissue. This virus is phloem limited and special care should be used during extraction.

Please see the user guide for full details.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)
- [Tech Note](#)



Reagent Set for Barley yellow dwarf virus-PAV (BYDV-PAV), 96 testwells, SRA 27500/0096

Reagent Set for Barley yellow dwarf virus-PAV (BYDV-PAV), 500 testwells, SRA 27500/0500

Reagent Set for Barley yellow dwarf virus-PAV (BYDV-PAV), 1000 testwells, SRA 27500/1000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com



Life Technologies™

Your Molecular & Cell Technology Partner



agdia®

Reagent Set for Barley yellow dwarf virus-PAV (BYDV-PAV), 5000 testwells, SRA 27500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com