

ELISA Reagent Set for Grapevine fanleaf virus (GFLV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Grapevine fanleaf virus*(GFLV) detection in plant foliage.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- Tween 20
- General extract buffer 3, GEB3
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)
- [Tech Note](#)



Reagent Set for Grapevine fanleaf virus (GFLV), 96 testwells,SRA 27703/0096

Reagent Set for Grapevine fanleaf virus (GFLV), 500 testwells,SRA 27703/0500

Reagent Set for Grapevine fanleaf virus (GFLV), 1000 testwells,SRA 27703/1000

Reagent Set for Grapevine fanleaf virus (GFLV), 5000 testwell,SRA 27703/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com