

ELISA Reagent Set for Tomato aspermy virus (TAV)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This test is a qualitative serological assay for detection of *Tomato aspermy virus*(TAV) in plant foliage.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent
- Tween 20
- General extract buffer 3 (GEB3), powder
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Tomato aspermy virus (TAV), 96 testwells, SRA 37001/0096

Reagent Set for Tomato aspermy virus (TAV), 500 testwells, SRA 37001/0500

Reagent Set for Tomato aspermy virus (TAV), 1000 testwells, SRA 37001/1000

Reagent Set for Tomato aspermy virus (TAV), 5000 testwells, SRA 37001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com