

ELISA Reagent Set for Potato virus S (PVS)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal


Intended Use:

This test can be used to detect all strains of *Potato virus S* (PVS), so far reported, in potato foliage and tuber sprouts.

A weak reaction with some PVM isolates has been observed.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Potato virus S (PVS), 96 testwells, SRA 40000/0096

Reagent Set for Potato virus S (PVS), 500 testwells, SRA 40000/0500

Reagent Set for Potato virus S (PVS), 1000 testwells, SRA 40000/1000

Reagent Set for Potato virus S (PVS), 5000 testwells, SRA 40000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com