

ELISA Reagent Set for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The ELISA is a reliable tool for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm), the causal agent of bacterial canker in tomato. The test can be used with bacterial cultures and ground leaf or stem tissue.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm), 96 testwells, SRA 44001/0096

Reagent Set for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm), 500 testwells, SRA 44001/0500

Reagent Set for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm), 1000 testwells, SRA 44001/1000

Reagent Set for *Clavibacter michiganensis* subsp. *michiganensis* (Cmm), 1000 testwells, SRA 44001/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com