

ELISA Reagent Set for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt)

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Synonym: *Corynebacterium michiganense* subsp. *tessellarium*

Intended Use:

The ELISA is a reliable tool for detecting *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), the causal agent of bacterial mosaic of wheat. The test can be used with bacterial cultures, wheat seed, or plant tissue.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), 96 testwells, SRA 44100/0096

Reagent Set for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), 500 testwells, SRA 44100/0500

Reagent Set for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), 1000 testwells, SRA 44100/1000

Reagent Set for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), 5000 testwells, SRA 44100/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com