

## ELISA Reagent Set for DsMV

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

---

This ELISA test is a qualitative serological assay for detection of *Dasheen mosaic virus* (DsMV) in cultivated species of the family Araceae.

Severe symptoms can result from DsMV infection in some plants. Mosaic and leaf distortion are common symptoms in crops such as dieffenbachia. *Dasheen mosaic virus* infection can also result in yield losses of up to 60 percent in such crops.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Dasheen mosaic virus (DsMV), 96 testwells, SRA 48000/0096

Reagent Set for Dasheen mosaic virus (DsMV), 500 testwells, SRA 48000/0500

Reagent Set for Dasheen mosaic virus (DsMV), 1000 testwells, SRA 48000/1000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)



**Life Technologies™**

*Your Molecular & Cell Technology Partner*



Reagent Set for Dasheen mosaic virus (DsMV), 5000 testwells, SRA 48000/5000

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)