

ELISA Reagent Set for PVA

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This test is used for detection of all strains of *Potato virus A* (PVA), so far reported, in potato foliage and sprouts.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Potato virus A (PVA), 96 testwells, SRA 60000/0096

Reagent Set for Potato virus A (PVA), 500 testwells, SRA 60000/0500

Reagent Set for Potato virus A (PVA), 1000 testwells, SRA 60000/1000

Reagent Set for Potato virus A (PVA), 5000 testwells, SRA 60000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com