

ELISA Reagent Set for CNFV

Test Label: Alkaline Phosphatase

Test Format: Antigen coated plate (ACP) ELISA

Capture Reagent: Polystyrene Plate

Detection Reagent: Monoclonal with Anti-animal conjugate

Intended Use:

This ELISA is a qualitative serological assay for detection of *Carnation necrotic fleck virus* (CNFV) in carnation and dianthus.

Includes:

- Detection antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Carnation necrotic fleck virus (CNFV), 96 testwells, SRA 68500/0096

Reagent Set for Carnation necrotic fleck virus (CNFV), 500 testwells, SRA 68500/0500

Reagent Set for Carnation necrotic fleck virus (CNFV), 1000 testwells, SRA 68500/1000

Reagent Set for Carnation necrotic fleck virus (CNFV), 5000 testwells, SRA 68500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com