

ELISA Reagent Set for AV-2

Test Label: Alkaline Phosphatase

Test Format: Antigen coated plate (ACP) ELISA

Capture Reagent: Polystyrene Plate

Detection Reagent: Polyclonal

Intended Use:

This test is used to detect AV2 in asparagus.

Includes:

- Detection antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Asparagus virus 2 (AV2), 96 testwells,SRA 71000/0096

Reagent Set for Asparagus virus 2 (AV2), 500 testwells,SRA 71000/0500

Reagent Set for Asparagus virus 2 (AV2), 1000 testwells,SRA 71000/1000

Reagent Set for Asparagus virus 2 (AV2), 5000 testwells,SRA 71000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com