

ELISA Reagent Set for CymMV

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This ELISA test is a qualitative serological assay for detection of *Cymbidium mosaic virus* (CymMV) in orchids. This ELISA test can detect CymMV in asymptomatic plants.

Cymbidium mosaic virus (CymMV) is one of the most prevalent viruses found in cultivated orchids. The virus is spread throughout orchid collections by way of mechanical transmission such as contaminated tools and pots. There is no known insect vector for this pathogen.

This test has been validated to detect CymMV isolates from Japan, Korea, and other parts of the world. No cross-reactions were observed when this assay was tested against other common pathogens of orchid (CMV, CymRSV, INSV, ORSV, Poty, TMV, ToRSV, TSWV) or related *Potexviruses* (PVX, HVX, PapMV).

If you are testing orchid, you may also be interested in:


Agdia's ELISA for detection of *Odontoglossum ringspot virus* (ORSV), another predominant pathogen in orchid cultivars (see related products for link).

Agdia's ImmunoStrip® for the detection of CymMV and ORSV: a simple and rapid testing device that can independently test for both pathogens. Perfect for use by hobbyists or greenhouse growers wishing to spot check their plants (see related products for link).

Testing services for orchid pathogens: Ask about our [orchid crop screen!](#) You can send your samples directly to our laboratory to be tested against an entire panel of pathogens known to infect orchid.

Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)
- [Tech Note](#)



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com

Reagent Set for Cymbidium mosaic virus (CymMV), 96 testwells, SRA 89501/0096

Reagent Set for Cymbidium mosaic virus (CymMV), 500 testwells, SRA 89501/0500

Reagent Set for Cymbidium mosaic virus (CymMV), 1000 testwells, SRA 89501/1000

Reagent Set for Cymbidium mosaic virus (CymMV), 5000 testwells, SRA 89501/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com