

ELISA Reagent Set for PnMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal


Intended Use:

This ELISA test is a qualitative serological assay for detection of *Poinsettia mosaic virus* (PnMV) in poinsettia plants.

PnMV is found in most of the branched, commercial poinsettia cultivars.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Poinsettia mosaic virus (PnMV), 96 testwells, SRA 90700/0096

Reagent Set for Poinsettia mosaic virus (PnMV), 500 testwells, SRA 90700/0500

Reagent Set for Poinsettia mosaic virus (PnMV), 1000 testwells, SRA 90700/1000

Reagent Set for Poinsettia mosaic virus (PnMV), 5000 testwells, SRA 90700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com