

ELISA Reagent Set for PDV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Prune dwarf virus* (PDV) in stone fruit crops.

ELISA has proven to be an excellent diagnostic method for this pathogen as it reduces the time and labor involved with older techniques, such as grafting.

Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- RUB6 enzyme conjugate diluent, 1X ready to use
- Tween 20
- General extract buffer 3, GEB3
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



Reagent Set for Prune dwarf virus (PDV), 96 testwells, SRA 98700/0096

Reagent Set for Prune dwarf virus (PDV), 500 testwells, SRA 98700/0500

Reagent Set for Prune dwarf virus (PDV), 1000 testwells, SRA 98700/1000

Reagent Set for Prune dwarf virus (PDV), 5000 testwells, SRA 98700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com