

ELISA Reagent Set for CaMV

Test Label: Peroxidase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The CaMV ELISA is a reliable method of detection of *Cauliflower mosaic virus* in foliage.

Includes:

- Capture antibody
- Peroxidase enzyme conjugate
- TMB peroxidase substrate solution, 1X
- MRS component, 5X
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)



ELISA Reagent Set, Cauliflower mosaic virus (CaMV),SRP 45600/0096

ELISA Reagent Set, Cauliflower mosaic virus (CaMV),SRP 45600/0500

ELISA Reagent Set, Cauliflower mosaic virus (CaMV),SRP 45600/1000

ELISA Reagent Set, Cauliflower mosaic virus (CaMV),SRP 45600/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com