

Test Tube Rack, 40 Places, Stainless Steel

Product ID T171

Description

Stainless steel wire test tube rack with 40-positions (4 x 10 rows). Holds 25 mm test tubes (products C930 and C935).

Product Information

Autoclavability	Up to 121 °C
Use With	C930, C935
Material	Stainless Steel
Places	40
Product Number	T171
Length	300 mm
Width	125 mm
Height	55 mm

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com