

Test Tube Rack, PhytoTech® Brand, 40-places, Wire, 25 mm

Product ID T8080

Description

A PhytoTech® brand product.

This is an epoxy coated wire test tube rack for 25 mm diameter test tubes. It has 40 places

Product Information

Autoclavability	Up to 121 °C
Material	Epoxy Coated Wire
Places	40
Length	285 mm
Width	122 mm
Height	88 mm

Shipping and Storage

Tariff Code	3926.90.9910
-------------	--------------

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com