

meta-Topolin


Product ID: **T841**

Introduction

Description: *meta*-Topolin is an adenine-based cytokinin that generally causes less root inhibition than B800, [6-(3-Hydroxybenzylamino)Purine

Product Information

Solubility	KOH
Physical Form	Solid
Storage Temp.	-20 to 0 °C
Other Notes	Plant Tissue Culture Tested
UPC / SKU	T841
CAS NUMBER	75737-38-1
Formula Weight	241.25
Formula	C ₁₂ H ₁₀ N ₅ OH
Synonyms	6-(3-Hydroxybenzylamino)purine



India Contact:



Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, opposite M2K Pitampura, Delhi-110034 (India)
Tel # +91-11-4220-8000; 4220-8111; 4220-8222 Fax# +91-11-4220-8444, Mobile# +91-98105-21400,
Toll Free # 1800-120-2434
Email - customerservice@lifetechindia.com | customerservice@atzlabs.com
ISO Certifications: [ISO 9001:2008 \(QMS\)](#) | [ISO 10002:2014](#)

Life Technologies, Lifetechindia.com, Lifetechindia.co.in , Lifetechnologiesindia.com , Lifetechnologiesindia.co.in, LifeTech, LifeTech (triangular arc), Auprep, Auprep.com, AuPreP Oligos Gelrose, Novascript, neofect**AMINE 3000**, genfect**AMINE 4000**, ATZ labs, ATZ Touch, ATZ Dew, ATZ Fabrich, ATZ MobiFriend, ATZLabs.com are the key trademarks/domains of Life Technologies India Private Limited and its associated entities.

Important: We have not granted licence or permission to any person or entity to use our above trade marks / domains
Any person / entity found to be using our marks shall be liable for civil / criminal action.

All Products shown on our website are for Laboratory research applications unless otherwise stated

© 2014. Life Technologies (India) Pvt. Ltd. All rights reserved.
design by www.nextconcept.in