

## Certificate of Analysis

### Human Transferrin (Tf) Protein, Holo

Cat. # TF25-N-BLK

Human Transferrin, Holo (Plasma)

**SIZE:** Bulk

**Lot#** TF2008-05

**FORM:** Soln

Lyophilized

Elemental iron is required for a variety of normal cellular functions and vital for proper growth and development. However, natural iron is quite insoluble and excess iron is harmful, since it can catalyze the formation of potentially damaging reactive oxygen species. The major pool of body iron (~85%; 40-50 mg/kg) is found in circulating hemoglobin and muscle myoglobin. Iron absorption occurs primarily in the intestine (duodenum) and inversely related to body iron reserve. Several proteins including **Ferritin, transferrin (Tf), transferrin receptors (TfRs)**, and **iron regulatory proteins (IRPs)** etc play a key role in iron metabolism.

**Transferrin** (Tf, human chromosome 3, 679 aa), a serum glycoprotein of ~80 kDa and synthesized in the liver, is the primary protein of inter-organ transport of nonheme iron. Tf can bind two iron atoms. Tf binds to membrane **Transferrin receptors (TfRs)** and taken up by endocytosis. Iron is released from Tf, within acidic endosomes, into the cytoplasm apparently through the action of DMT1. The apoTf-TfR complex is returned to the cell surface, where, apo-Tf dissociates from TfR at the extracellular pH. The classical TfR, now termed **TfR1**, is a homodimeric (95 kDa subunits) type II membrane glycoprotein that binds two molecules of Tf. Human TfR1 (human 760 aa; mouse 763 aa) has a cytoplasmic domain 1-67aa, 68-88 aa TM, and 89-760 aa as extracellular domains. A monomeric serum form or **soluble TfR1** (~80 kDa) also exists that lacks residues 1-100 aa. Recently, a second Tf receptor, **TfR2**, has been cloned and characterized. TfR2 shares 45% identity with TfR1, and 27% with PMSA. Several variants of Tf have been identified with varying iron binding ability.

#### Source of Antigen

Prepared from plasma shown to be non-reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA required tests. US origin.

**Purity:** Purified human plasma Tf (>97% purity by SDS-PAGE)

**Form:** Each gram lyophilized from 7.0 ml 20 mM phosphate buffer, pH 7.4, 150 mM NaCl

**Protein Concentration:** 143 mg/ml determined by Biuret

**Iron content:** >=1.0ug/mg protein; Total Iron Binding Capacity: 1.2 to 1.6 ug/mg

**Storage:** -20oC or below

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for powder.

**General References:** Bowman, B. H. et al (1988) Adv. Genet. 25: 1-38; Evans, R. W. et al (1982) Biochem. J. 201: 19-26; MacGillivray, R. T. A et al (1982) PNAS 79: 2504-2508; Park, I. et al (1985) PNAS 82, 3149; Uzan, G. et al (1984) BBRC 119, 273; Yang, F. et al (1984) PNAS 81, 2752-2756; Nelson N (1999) EMBO J. 18, 4361; Cairo G I (2000) Biochem. J. 352, 241-250

THIS PRODUCT IS A LABORATORY REAGENT.  
IT IS NOT TO BE ADMINISTERED TO HUMANS  
OR USED FOR ANY DRUG PURPOSE.  
FOR RESEARCH USE ONLY

#### Related material available from ADI

Antibodies NRAMP1/2, MTP1, Transferrin, and TfRs receptors (TfR1 and TfR2), Ferritin, H and L-chain, HOX1-3, HFE, Dcytb, IRP1 and IRP2, Frataxin; Human and Mouse Transferrin ELISA Kits

TF25-N-BLK

90713A

#### India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)