

**Product Specification Sheet**

**Troponin-I Sheep Skeletal Muscle**

**Cat. #** TNIM45-N-25      Troponin-I Sheep Skeletal Muscle, pure protein      **SIZE:** 25 ug

Troponins are protein components of striated muscle. There are three different troponins: troponin C, troponin T and troponin I. Troponin T (cTnT) and troponin I (cTnI) are released only following cardiac damage. They are present for, and remain elevated, a long time. Unlike CK and CK-MB, cTnT and cTnI are released for much longer with cTnI detectable in the blood for up to 5 days and cTnT for 7-10 days following MI. This allows an MI to be detected if the patient presents late.

Troponin T and I are very sensitive. There is always a low level release of CK and CK-MB from skeletal muscle at a low level all the time so there is always a background value. This is not the case for the cardiac structural proteins such as cTnT and cTnI and therefore, they are very sensitive. Studies have revealed that about one third of patients admitted with unstable angina, in which MI was apparently excluded by CK and CK-MB measurement, have raised levels of cTnT and cTnI.

Troponin C, skeletal muscle is a protein that in humans is encoded by the TNNC2 gene. Troponin (Tn), a key protein complex in the regulation of striated muscle contraction, is composed of 3 subunits. The Tn-I subunit inhibits actomyosin ATPase, the Tn-T subunit binds tropomyosin and Tn-C, while the Tn-C subunit binds calcium and overcomes the inhibitory action of the troponin complex on actin filaments. The protein encoded by this gene is the Tn-C subunit.

cTnI is a part of the troponin complex. It binds to actin in thin myofilaments to hold the actin-tropomyosin complex in place. Because of it myosin cannot bind actin in relaxed muscle. When calcium binds to the Troponin C it causes conformational changes which lead to dislocation of troponin I and finally tropomyosin leaves the binding site for myosin on actin leading to contraction of muscle. The letter I is given due to its inhibitory character.

**Source of Antigen and Antibodies**

Monkey Troponin I was purified from African green Monkey skeletal muscle using proprietary techniques (>95%, ~21 kda). Purified protein (cat # TNIM45-N-25) is formulated in 8M Urea, 50 mM Tris.(pH 8.0), 2 mM EDTA, and 5mM beta-mercaptoethanol (see lot sp conc on the vial).

This preparation can be used for ELISA or diluted in other appropriate buffers.

Store frozen in suitable aliquots at -20 °C or lower.

**Purified protein**

- 25 ug/vial
- solution
- lyophilized powder

**Reconstitute** powder in appropriate buffer in at least 100 ug/ml.

**Storage**

**Short-term:** unopened, undiluted vials for less than a week at 4oC.

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for lyophilized items.

**General References:** (1).

\*This product is for in vitro research use only.

**Related items**

Catalog#	ProdDescription
TMY16-N-25	Tropomyosin, Dog Cardiac, purified protein
TMY17-N-25	Tropomyosin, Human Cardiac, purified protein
TNCC75-N-25	Troponin-C, Human Cardiac, pure protein
TNCC76-N-25	Troponin-C, Recombinant, Human Cardiac, pure protein
TNIC15-A	Anti-Human Troponin-I Cardiac, (peptide 3), Monoclonal Anti-Human Troponin-I, Cardiac,
TNIC16-M	Troponin-I, Recombinant, Human Cardiac, pure protein
TNIC25-N-25	Troponin-I, Bovine Cardiac, pure protein
TNIC26-N-25	Troponin-I, Mouse Cardiac, pure protein
TNIC27-N-25	Troponin-I Pig Cardiac, pure protein
TNIC28-N-25	Troponin-I Pig Cardiac, pure protein
TNIC29-N-25	Troponin-ITC Complex Pig Cardiac , pure protein
TNIC31-N-25	Troponin-I Rabbit Cardiac, pure protein
TNIC32-N-25	Troponin-I Rat Cardiac, pure protein
TNIC33-N-25	Troponin-I, Dog Cardiac, pure protein
TNIM49-N-25	Troponin-C Rabbit Skeletal Muscle, pure protein
TNIM47-N-25	Troponin-I Rat Skeletal Muscle, pure protein
TNIM51-N-25	Troponin-I, Chicken Skeletal Muscle, pure protein
TNIM53-N-25	Troponin-I Human Skeletal Muscle, pure protein
TNTC55-A	Monoclonal Anti-Human Troponin-T, Cardiac,
TNTC65-N-25	Troponin-T, Human Cardiac, pure protein

TNIM45-N-25

110124A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)